# GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy

Italia: tel. 199 400 401 - fax 199 400 403

Export: tel. +39 02 953854209/221/225 fax +39 02 95380056

[gima@gimaitaly.com](mailto:gima@gimaitaly.com) - [export@gimaitaly.com](mailto:export@gimaitaly.com)

[www.gimaitaly.com](http://www.gimaitaly.com)

---

**SATURIMETRO PEDIATRICO OXY-PED**  
**OXY-PED PEDIATRIC FINGERTIP OXIMETER**  
**OXYMÈTRE PÉDIATRIQUE OXY-PED**  
**PULSOXIMETER FÜR KINDER OXY-PED**  
**SATURÓMETRO PEDIÁTRICO OXY-PED**  
**MEDIDOR DE SATURAÇÃO PEDIÁTRICO OXY-PED**  
**ΠΑΙΔΙΑΤΡΙΚΟ ΚΟΡΕΣΤΟΜΕΤΡΟ OXY-PED**

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## Instructions to User

Read these instructions carefully before using this equipment. These instructions describe the operating procedures to be followed strictly. Failure to follow these instructions can cause measuring abnormality, equipment damage and personal injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to user's negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

Please peruse the relative content about the clinical restrictions and caution.

### 1. Safety

#### 1.1 Instructions for Safe Operations

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance. It is recommended that the device should be inspected once a week at least. When there is obvious damage, stop using the oximeter.
- Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
- The oximeter cannot be used together with the devices not specified in User's Manual.

#### 1.2 Warnings

- *Explosive hazard—DO NOT use the oximeter in environment with inflammable gas such as some ignitable anesthetic agents.*
- *DO NOT use the oximeter while the testee measured by MRI and CT*
- *To dispose the device, the local law must be followed.*

### 1.3 Attentions

- *Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.*
- *If the oximeter gets wet, please stop operating it. When it is carried from cold environment to warm and humid environment, please do not use it immediately.*
- *DO NOT press the keys on front panel with sharp materials.*
- *High temperature or high pressure steam disinfection to the oximeter is not permitted. Refer to User's Manual for instructions of cleaning and disinfection.*
- *DO NOT have the oximeter immersed in liquid. When it needs cleaning, please wipe its surface with disinfect solution by soft material. Do not spray any liquid on the device directly.*

## 2. Overview

The pulse oxygen saturation is the percentage of  $\text{HbO}_2$  in the total Hb of the blood, so-called the  $\text{O}_2$  concentration in the blood. It is an important bio-parameter to the respiration. Many of the respiration disease will cause hypoxemia, even endanger the patient's health. As a result, monitoring the  $\text{SpO}_2$  is indispensable in the clinical rescuing. The traditional method to measure  $\text{SpO}_2$  is to analyze the sample of the patient's blood, so can get the partial pressure of oxygen and calculate the  $\text{SpO}_2$  by use the blood-gas analyzer. This method can not be used to monitor continuously. For the purpose of measuring the  $\text{SpO}_2$  more easily and accurately, GIMA developed the Fingertip Oximeter. The device can measure the pulse rate and blood perfusion index simultaneously.

The Fingertip Oximeter is compact, convenient to use and carry and with low power consumption. You just need to put the fingertip into the sensor of the device, the  $\text{SpO}_2$  value will appear on the screen immediately.

### 2.1 Features

- Small in size and lightweight;
- Color OLED, various display modes, display directions adjustable;

- The device can accurately measure SpO<sub>2</sub> value, pulse rate value and perfusion index;
- The device will automatically start measuring after putting finger into rubber cushions;
- The device will power off automatically without signal for about 8 seconds;
- Audible & visual alarm function;
- Low voltage indication;

## 2.2 Major Applications and Scope

This device is applicable to home, hospital (including internal medicine, surgery, anesthesia, pediatrics, emergency room etc.), oxygen bar, the community medical center, alpine area and it also can be used before or after sports, and the like. The Fingertip Oximeter can detect SpO<sub>2</sub>, pulse rate and blood perfusion index.

*This device is not appropriate to be used for continuous monitoring.*

## 2.3 Environment Requirements

Operating Temperature: 5°C ~40°C

Operating Humidity: 30%~80%

Atmospheric pressure: 70kPa~106kPa

## 2.4 Principle of Measurement

The measurement of pulse oximeter is that it uses a multi-functional oxyhemoglobinometer to transmit some narrow spectrum light bands through blood samples and to measure attenuation of spectrum with different wavelengths according to the characteristic that RHb, O<sub>2</sub>Hb, Met Hb and COHb absorb the light of different wavelength, thereby determining O<sub>2</sub>Hb saturation of different fractions. O<sub>2</sub>Hb saturation is called "fractional" O<sub>2</sub>Hb saturation.

Fractional O<sub>2</sub>Hb saturation =  $[O_2Hb / (RHb + O_2Hb + Met Hb + COHb)] \times 100$

Oppositely, pulse oxygen oximeter measures functional O<sub>2</sub>Hb saturation:

$$\text{Functional O}_2\text{Hb saturation} = [\text{O}_2\text{Hb} / (\text{RHb} + \text{O}_2\text{Hb})] \times 100$$

Present SpO<sub>2</sub> oximeter transmits light of two wavelengths only, red light (wavelength 660 nm) and infrared (wavelength 940 nm), to differentiate HbO<sub>2</sub> from HbR. One side of the sensor contains two LEDs, and the other side contains a photoelectric detector. SpO<sub>2</sub> oximeter measures HbO<sub>2</sub> saturation in the blood by the light plethysmograph when the pulse beats. The result is quite precise when HbO<sub>2</sub> saturation is over 70% ~ 95%.

### 2.5 Caution

A. The uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.

B. For the individual patients, there should be a more prudent inspecting in the placing process. The device can not be clipped on the edema and tender tissue.

C. The light (the infrared is invisible) emitted from the device is harmful to the eyes, so the user and the maintenance man, can not stare at the light.

D. Testee can not use enamel or other makeup on the finger.

E. Testee's fingernail can not be too long.

F. The finger should be placed properly, or else it may cause inaccurate measurement.

G. Do not shake the finger and keep at ease during using.

H. Do not put wet finger directly into oximeter.

I. The SpO<sub>2</sub> sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.

L. Do not fix the SpO<sub>2</sub> sensor with adhesive or else it may result in venous pulsation and inaccurate measure of SpO<sub>2</sub>.

M. Make sure the optical path is free from any optical obstacles like rubberized fabric.

N. Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight and etc.

O. Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.

### 3. Technical Specifications

**A. Display mode:** OLED Display

**B. Power supply requirement:**

Rechargeable lithium battery: one coin cell

Battery voltage: 3.6VDC

**C. Operating current:**  $\leq 50\text{mA}$

**D. SpO<sub>2</sub> Parameter Specifications:**

Measuring range: 35% - 99%

Accuracy:  $\pm 2\%$  (during 75% - 99%)

Accuracy:  $\pm 3\%$  (during 50% - 74%)

**E. Pulse Rate Parameter Specifications:**

Measuring range: 30bpm-240bpm

Accuracy:  $\pm 2\text{bpm}$  or  $\pm 2\%$  (which ever is greater)

**F. Blood Perfusion Parameter Specifications:**

Measuring range: 0%~20%

**G. Resistance to interference of surrounding light:**

The difference between the value measured in the condition of indoor natural light and that of darkroom is less than  $\pm 1\%$ .

**H. Resistance to interference of man-made light:**

Values of SpO<sub>2</sub> and Pulse Rate can be accurately measured by pulse oxygen simulator.

**I. Dimensions:** 56 mm (L) x 47 mm (W) x 32 mm (H) **Net Weight:** 45g (including batteries)

**J. Classification:**

The type of protection against electric shock: Internally powered equipment.

The degree of protection against electric shock: Type BF applied part.

The degree of protection against harmful ingress of liquids: Ordinary equipment without protection against ingress of water.

Electro-Magnetic Compatibility: Group I, Class B


#### 4. Accessories

- A. A hanging cord
- B. Two batteries
- C. A charger
- D. A power adapter
- E. A charger cable
- F. A pouch
- G. A User Manual
- H. Quality Certificate

#### 5. Installation

##### 5.1. Front view (Figure 1)

##### 5.2. Battery

1. Insert the coin cell into the battery compartment properly in the right direction (as shown in Figure 2).


Figure 2


2. Replace the cover.

- Please take care when you insert the coin cell, improper insertion may make the oximeter not work.

## 6. Operating Guide

1. Open the clip as shown in Figure 3.

Figure 3

2. Put finger onto the rubber cushion of the clip (make sure the finger is in the right position), and then clip the finger.

3. The device will power on automatically in 2 seconds, and start to display software version number.

4. Next enter into data display screen (as shown in Figure 4). User can read the values and view the waveform from display screen.


Figure 4 PR display


Figure 5 PI% display

Screen description:

“%SpO<sub>2</sub>”: SpO<sub>2</sub> symbol; “99”: SpO<sub>2</sub> value;

“PR”: Pulse rate icon; “65”: Pulse rate value;

“♥”: Pulse beat symbol;

“T”: Pulse intensity histogram.

“PI”: Perfusion index icon; “1.4”: Perfusion index value;

5. At this time, longtime press Display Key (about 2 seconds), display screen will be shown as Figure 5. The PR is replaced by PI%. To return to PR display, just longtime press Display key again. If there is no operation in 20 seconds, it will return to PR display as well.

6. Change display direction.

Four directions display alternately. Press Display Key once at a time, display screen will be flipped 90° each time in a cyclic manner. But there is a little difference, the display facing left or right will not show plethysmogram.

7. Alarm Indicator

When measuring, if SpO<sub>2</sub> value and pulse rate value exceeds the preset alarm limits, the device will alarm automatically and the value on the screen exceeding limit will blink; at this time press Display Key to shut down the alarm. Exceeding pulse rate alarm limit: sound twice as an interval.

Preset alarm range:


SpO<sub>2</sub> alarm: Low limit: 90%

Pulse Rate alarm: High limit 120bpm - Low limit: 50bpm


## 7. Instructions for Charging the Cell

The rechargeable coin cell is used to supply power for the oximeter. When it is run out, you need to use the provided charger to charge it, details is as follows:

1. Insert the cell into charger; pay attention to its polarity (as shown in the following figure).

2. Connect the charger and power adapter by charger cable in the following way.

3. The power adapter should be plugged into the AC mains with 100~250VAC, 50/60Hz.

4. Status of charging indicator

Green: full

Orange: no cell

Red: charging

Note:

- The charger will stop charging automatically when the coin cell is full (the indicator becomes green). Moreover, it has over heat protection function, when the coin cell temperature is higher than 45°C, it will also stop charging automatically.
- The charger can be connected to the USB port of a computer by the charger cable for charging.

## 8. Repair and Maintenance

- Please charge the battery when the low-voltage indicator lightens.
- Please clean the surface of the device before using. Wipe the device with alcohol first, and then let it dry in air or wipe it dry.
- Please take out the battery if the oximeter will not be used for a long time.
- The best storage environment of the device is -20°C to 55°C ambient temperature and less than 95% relative humidity.

*High-pressure sterilization cannot be used on the device.*

*Do not immerse the device in liquid.*

*It is recommended that the device should be kept in a dry environment. Humidity may reduce the useful life of the device, or even damage it.*

### 8.1 Battery Maintenance

Keep the both sides of coin cell clean.

- Low temperature may decrease the performance of coin cell, and low battery indicator may appear early. In such case, please put coin cell into pocket for warm before use, thus bring it back to normal condition.
- Do not let any conductive metal (such as tweezers) contact both sides of coin cell simultaneously to avoid short circuit.
- Charge the coin cell for 8~10 hours each time; ambient temperature should be 5°C~40°C.
- If the coin cell is full after charging, but its performance decreases apparently, it means the coin cell is exhausted, please change a new one.

## 9. Troubleshooting

Trouble	Possible reason	Solution
<b>The SpO<sub>2</sub> and Pulse Rate display instable</b>	<ol style="list-style-type: none"><li>1. The finger is not placed inside enough.</li><li>2. The finger is shaking or the patient is moving.</li></ol>	<ol style="list-style-type: none"><li>1. Place the finger properly and try again.</li><li>2. Let the patient keep calm.</li></ol>
<b>The device can not turn on</b>	<ol style="list-style-type: none"><li>1. The coin cell is drained or almost drained.</li><li>2. The coin cell is not inserted properly.</li><li>3. The device's malfunction.</li></ol>	<ol style="list-style-type: none"><li>1. Charge the coin cell.</li><li>2. Reinstall coin cell.</li><li>3. Please contact the local service center.</li></ol>
<b>The indicator light is off suddenly</b>	<ol style="list-style-type: none"><li>1. The device will power off automatically when it gets no signal for 8 secs.</li><li>2. The coin cell is almost drained.</li></ol>	<ol style="list-style-type: none"><li>1. Normal.</li><li>2. Charge the coin cell.</li></ol>

**Disposal:** *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.*

### **GIMA WARRANTY CONDITIONS**

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons.

Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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