GIMA

**ASPIRATORE VEGA
ASPIRATOR VEGA
ASPIRATEUR VEGA
ABSAUGER VEGA
ASPIRADOR VEGA**

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VEGA it's a device working 230V ~ / 50 Hz network electricity, to be used for the nasal aspiration, oral aspiration, tracheal aspiration of the body liquids (mucus, catarrh or blood) in the adult or in the children. Easily portable equipment designed for not continuous use 30 min. ON / 30 min. OFF.
Thanks to this characteristics and to the rating that it has, this product is particularly suitable for hospital use, on the tracheotomized patients, minor surgical applications and post-operative therapy at home.
Made of highly heat resistant, electrically insulated plastic material in conformity with the latest European safety standard, the product is supplied with a complete polycarbonate autoclavable jar with overflow valve and it is equipped with aspiration regulator and vacuum indicator located on the front panel.

GENERAL WARNINGS

READ INSTRUCTION MANUAL CAREFULLY BEFORE USE

ONLY HIGHLY QUALIFIED STAFF USE RESERVED

**THE INSTRUMENT MUST NOT DISASSEMBLED
FOR TECHNICAL SERVICE ALWAYS CONTACT GIMA**

IMPORTANT SAFETY RULES

1. When opening packing verifies the device state: always check the presence of damages on the plastic surface, this can expose under tension internal parts. Also check the main cable integrity.
In case of damage do not connect the plug to the socket. Check integrity before every use.
2. Before connecting the instrument, always check that the electrical requirements on manufacturer label and the plug used comply with the electrical network used.
3. Keep special attention to:
 - Use only original accessories;
 - The device can be used only with the bacteriological filter;
 - Place instrument on stable and flat surfaces and to avoid the blockage of the air inlets on sides;
 - Do not use in presence of air inflammable anaesthetic, oxygen or nitrous oxide;
 - Do not touch the instrument with we hands, and avoid contact with liquids;
 - Keep off the reach of children or not capable people without supervision;
 - If not used do not leave the instrument connected to network;
 - Do not pull cable to remove plug from socket; disconnect hold plug with fingers;
 - Use and keep the instrument in a safe environment, protected from bad weather condition and keep off excessive heating;
4. **The device must be used for the use it was designed.**
A different use must be considered improper and dangerous. Manufacture cannot be considered responsible for damages due to improper use or connection to an electrical system not complying with current regulations. **Do not use for anything else but the use defined by manufacturer;**
5. Instrument and accessories discharging must be done following current law regulations in every country of use.
6. For repairs, exclusively contact technical service and request the use of original spare parts.
Failure to comply with the above can jeopardise the safety of the device;
7. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents;

IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:

At the end of its working life, the product must not be disposed of as urban waste.
It must be taken a special local authority differentiated waste collection or to a dealer providing this service.
Disposing of a household appliance separately avoids possible negative consequences for the environment and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant saving in energy and resources. As a reminder of the need to dispose of separately the product is marked with a crossed-out wheeled dustbin.

TECHNICAL CHARACTERISTICS

Model	VEGA	
Typology (MDD 93/42/EEC)	Medical Device Class II a	
Classification UNI EN ISO 10079-1	HIGH VACUUM / LOW FLOW	
Main Voltage	230 V ~ / 50 Hz	110 V ~ / 60 Hz (no CE 0123)
Power consumption	184 VA	105 VA
Fuse	F 1 x 1.6A 250 V	F 1 x 4 A 250 V
Maximum suction aspiration (without jar)	-75kPa (- 0.75 bar) Adjustable from -75kPa (-0.75bar) to -10kPa (-0.10 bar)	
Maximum flow (without jar)	16 l/min	
Weight	2.5 Kg	
Dimension	350 x 210 x 180 mm	
Functioning (to 35°C and 110 % operating voltage)	30 min ON / 30 min OFF	
Working condition	Room temperature: 5 ± 35 °C Room humidity percentage: 30 ± 75 % RH Altitude: 0 ± 2000m s.l.m	
Conservation condition	Room temperature: - 40 ± 70 °C Room humidity percentage: 10 ± 100% RH	

SYMBOLS

	Class 2 isolation equipment
CE 0123	CE marking in conformity with EC directive 93/42/CEE Manufacturer : CA-MI snc Via Ugo La Malfa n°31 43010 Pilaastro (PR) Italy
	Warning, consult the instruction manual
	Type B equipment
	Fuse
	To preserve in place coolness and dry land
	Conservation temperature: - 40 ± 70° C
~	Alternate current
Hz	Mains frequency
I	ON
0	OFF

Guidance and manufacturer's declaration – electromagnetic Emissions (as request by regulation EN 60601-1-2:2001)		
The surgical aspirator VEGA is intended for use in the electromagnetic environment specified below. The customer or the user of the surgical aspirator VEGA should assure that it's used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
Irradiated / Conducted emissions CISPR11	Group 1	The surgical aspirator VEGA only used RF energy only for its internal functioning. Therefore, its RF emissions are very low and are not cause interference in proximity of any Electronic appliances. The surgical aspirator VEGA can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes.
Irradiated / Conducted emissions CISPR11	Class [B]	
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic Immunity (as request by regulation EN 60601-1-2:2001)		
The surgical aspirator VEGA is intended for use in the electromagnetic environment specified below. The customer or the user of the surgical aspirator VEGA should assure that it's used in such an environment.		
Immunity Test	Compliance	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC/EN 61000-4-2	± 6kV on contact ± 8kV in air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	± 2kV power supply	Mains power quality should be that of a typical commercial environment or hospital.
Surge IEC/EN 61000-4-5	± 1kV differential mode	Mains power quality should be that of a typical commercial environment or hospital.
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	5%U _T for 0.5 cycle 40%U _T for 05 cycle 70%U _T for 25 cycle <5%U _T for 5 sec	Mains power quality should be that of a typical commercial environment or hospital. If the user of the surgical aspirator VEGA request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Conducted Immunity IEC/EN 61000-4-6	3Vrms 150kHz to 80MHz (for appliances that aren't life - supporting)	-
Irradiated Conducted IEC/EN 61000-4-3	3V/m 80MHz to 2.5 GHz (for appliances that aren't life – equipment)	-

Nota U_T is the value of the power supply voltage

ACCESSORIES SUPPLIES

DESCRIPTION
COMPLETE ASPIRATION JAR 1000cc
CONICAL FITTING
TUBES SET 6mm x 10mm (TRASPARENT SILICON)
ASPIRATION PROBE CH 20
ANTIBACTERIAL FILTER

The filter is produced with (PTFE) hydrophobic material witch prevents fluids entering the pneumatic circuit. When the filter is wet, it's not possible to use the unit therefore the filter should be changed immediately. In case of possible contamination or discolouration, change the filter immediately. Don't use the suction unit without the protection filter fitted. If the suction unit is used in an emergency or in a patient where the risk of contamination is not know the filter must be changed after each use. The suction probe must be changed at the end of each application. Available under request with different versions with two complete jar 2000cc.

CLEANING OF ACCESSORIES

Wearing disposable latex gloves, proceed to clean the accessories as follows:

1. Grip the jar with your hand and turn the cap in a counter-clockwise direction.
2. Remove the red cage covering the float located in the jar cap (conical part with 3 openings)
3. Once the cage is removed, the float with the red seal comes out of the jar cap.
4. **To disinfect the accessories, use either denatured alcohol or hypochlorate-based solutions, easily purchasable from chemist shops.**
5. Once cleaning has been completed, insert the float in the cage, with the seal facing towards the cage opening.
6. Insert the opening of the cage into the seat in the jar cap.
7. To facilitate the cleaning of the jar, fill it 1/3 full of ordinary water. The aspirated liquid, thus diluted, will then be easier to remove.

It is possible autoclave the accessories cover and jar: insert the parts in autoclave and to effect a cycle of sterilization with vapor to the temperature of 121°C (pressure relative 1 bar) having care to position turned upside-down the graduated jar (with turned fund toward the tall one).

After the sterilization and the cooling to the temperature the components' environment, to verify that this last don't result damaged; disassembled therefore the container for inhaled liquids following the inverse operations to the dismantlement.

The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 120°C. The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C. The device is ready for a new employment now.

During the sterilization operations the operator must use disposable gloves in latex, shirts and oportune mask not to enter contact with possible substances contaminating

DO NOT WASH, STERILIZE OR PUT IN AUTOCLAVE THE ANTIBACTERIAL FILTER

PERIODICAL MAINTENANCE CHECKS

The VEGA suction equipment does not need maintenance or lubrication. It is necessary to check functioning and instrument before every use. Unpack the instrument and **always check** integrity of plastic parts and feeding cable, they might have been damaged during previous use. Connect cable to electrical network and turn switch on. Close the aspirator outlet with your finger and with suction regulator in maximum vacuum position check that the vacuum indicators reaches -75kPa (-0.75 bar) minimum. Rotate the knob from right to left and check the aspiration regulating control. The vacuum indicator should go down -25kPa (-0.25 bar). Verify that loud noises are not present, these can indicate wrong functioning. A protection fuses (**F 1 x 1.6 A 250V or F 1x4A 250V for voltage 110V/60Hz**) reachable from exterior and situated in the plug protects the instrument. For fuses replacing, always the type and the range.

Type of fault	Cause	Remedy
1. No aspiration	Jar Cap badly screwed down	Unscrewed the cap, then rescrew it correctly
2. No aspiration	Lid seal not in its seat	Unscrew the cap and insert the seal properly in its seat
3. The float doesn't close	If the cap has been washed, ensure that the float is not partially detached	Insert the float into it's place
4. The float doesn't close	The float it's covered by dirty material	Unscrewed the cap, leave the and put in on autoclave
5. Low suction	Foam inside the jar	Fill the jar to 1/3 full of ordinary water
6. No aspiration due to flow leakage of mucus	Filter blocked	Replace filter
Faults 1 - 2 - 3 - 4 - 5 - 6	None of the remedies has achieved the desired results	Contact the seller or GIMA After-sales Assistance Service

If the overfill security system it's activated, don't proceede with the liquid aspiration. If the overfill security system doesn't work there are two cases:
 1° case – If the overfill security system doesn't work the aspiration will be stopped by the bacteriological filter who avoid the liquid penetration inside the device.
 2° case – If both the security system doesn't work, there is the possibility that liquid comes inside the device, in this case return the device to GIMA technical service.

**BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT GIMA TECHNICAL SERVICE.
 GIMA DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.**

INSTRUCTION FOR USE

- Connect the short silicon tube **1** ,with filter antibacterial **8**, to the suction connector **2**. The other tube, with one end connected to the filter must be connected with the other end to jar's lid connector **3** where has been fixed the red float. When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to avoid liquid penetration inside the device. The device must be used on a plan of horizontal operation.
- Connect the long silicon tube **4** to the other jar's lid connector **5**
- Connect the other end of the long silicon tube **4** to the probe plastic connector **6** then connect the suction probe to it.
- Connect the power cord to the device then connect the plug to the electrical mains supply.
- Push switch **7** on position **1** to start suction
- Once finished push switch on **0** position and unplug.
- Unscrew the jar's lid and fill the jar 1/3 full or ordinary water (this for an easy cleaning operations and an rapid reaching of the functionally vacuum) then rescrew the lid on the jar correctly.
- To extract the accessories and start with cleaning.

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